

IHE Integration Statement			
Vendor	Product Name	Version	Date
iCAD, Inc.	SecondLook® Digital / PowerLook AMP	7.2	9/25/15
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Mammography Image (MAMMO)	Evidence Creator	None	
Evidence Document (ED) Integration Profile	Evidence Creator	None	
Post-Processing Workflow (PWF) Integration Profile	Evidence Creator	None	
<u>Internet address for vendor's IHE information:</u> http://www.icadmed.com/assets/dtb084-ihe-integration-statement-sld.pdf			
Links to Standards Conformance Statements for the Implementation			
HL7	Not applicable		
DICOM	http://www.icadmed.com/assets/powerlook-amp-dicom-conformance-statement.pdf		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	